

### Recombinant Human RBP3/IRBP protein ab116892

画像数 1

#### 製品の詳細

製品名	Recombinant Human RBP3/IRBP protein
発現系	Wheat germ
アクセッション番号	<b><u>P10745</u></b>
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
配列	GTAEETYIMKRLGRALVIGEVTSGGCQPPQTYHVDDTNLYL TIPTARSV GASDGSSWEGVGVTPHVVPAAEEALARAKEMLQHNQLRVKRS PGLQDH
予測される分子量	36 kDa including tags
領域	1149 to 1246

#### 特性

Our **Abpromise guarantee** covers the use of **ab116892** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	ELISA Western blot SDS-PAGE
製品の状態	Liquid
備考	This product was previously labelled as RBP3.


#### 前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.3% Glutathione, 0.79% Tris HCl
------------	---

## 関連情報

機能	IRBP shuttles 11-cis and all trans retinoids between the retinol isomerase in the pigment epithelium and the visual pigments in the photoreceptor cells of the retina.
配列類似性	Belongs to the peptidase S41A family.
細胞内局在	Secreted > extracellular space > extracellular matrix > interphotoreceptor matrix. Interphotoreceptor matrix that permeates the space between the retina and the contiguous layer of pigment epithelium cells.

## 画像



12.5% SDS-PAGE showing ab116892 at approximately 36.41 kDa stained with Coomassie Blue

SDS-PAGE - Recombinant Human RBP3/IRBP protein (ab116892)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors