

Recombinant Human RCAS1 protein ab79243

画像数 1

製品の詳細

製品名	Recombinant Human RCAS1 protein
精製度	> 90 % SDS-PAGE. ab79243 is purified using conventional chromatography techniques.
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MGSSHHHHHH SSGLVPRGSH MRSGRGRKLS GDQITLPTTV DYSSVPKQTD VEEWTSWDED APTSVKIEGG NGNVATQQNS LEQLEPDYFK DMTPTIRKTQ KIVIKKREPL NFGIPDGSTG FSSRLAATQD LPFIHQSSSEL GDLDTWQENT NAWEEEEEDAA WQAEVLRQQ KLADREKRAA EQQRKKMEKE AQRLMKKEQN KIGVKLS

特性

Our **Abpromise guarantee** covers the use of **ab79243** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション SDS-PAGE

製品の状態 Liquid


前処理および保存

保存方法および安定性 Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
pH: 8.00
Constituents: 0.0154% DTT, 0.242% Tris, 0.0584% EDTA, 10% Glycerol (glycerin, glycerine)

関連情報

機能	May participate in suppression of cell proliferation and induces apoptotic cell death through activation of interleukin-1-beta converting enzyme (ICE)-like proteases.
組織特異性	Widely expressed. Expressed in ovary, testis, prostate, thymus, muscle and heart, but not in small intestine, colon, lymph nodes, or peripheral blood lymphocytes. The protein is not detected in any of the above organs.
ドメイン	The coiled coil domain is necessary for the homodimerization.
細胞内局在	Golgi apparatus membrane. According to PubMed:10426319, it also exists as a soluble form which has the same biological activities. The existence of such soluble form is however uncertain.

画像



15% SDS-PAGE showing ab79243 (3μg).

Note: Real molecular weight on SDS-PAGE is shifted up compared to predicted molecular weight.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors