

Recombinant Human RCN1/RCN protein ab98241

画像数 1

製品の詳細

製品名	Recombinant Human RCN1/RCN protein
精製度	>= 80 % SDS-PAGE. Purity > 80% by SDS-PAGE. ab98241 was purified using conventional chromatography techniques
発現系	Escherichia coli
アクセッション番号	<u>Q15293</u>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	<p>MRGSHHHHHHGMASMTGGQQMGRDLYDDDDKDRWGSEL EKPTVRKERVVR PDSELGERPPEDNQSFQYDHEAFLGKEDSKTFDQLTPDESKE RLGKIVDR IDNDGDGFVTTEELKTWIKRVQKRYIFDNVAKVWKDYDRDKD DKISWEEY KQATYGYLGNPAEFHDSSDHHTFKKMLPRDERRFKAADLNG DLTATREE FTAFLHPEEFEHMKEIVVLETLEDIDKNGDGFVDQDEYIADM FSHEENGP EPDWVLSEREQFNEFRDLNKDGKLDKDEIRHWILPQDYDHAQ AEARHLVY ESDKNKDEKLTKEEILENWNMFVGSQATNYGEDLTKNHDEL</p>
予測される分子量	40 kDa including tags
領域	30 to 331
タグ	His tag N-Terminus

特性

Our **Abpromise guarantee** covers the use of **ab98241** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション SDS-PAGE

	Mass Spectrometry
質量分析	MALDI-TOF
製品の状態	Liquid
備考	This product was previously labelled as RCN1


前処理および保存

保存方法および安定性	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. pH: 8.00 Constituents: 0.0154% DTT, 0.316% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.58% Sodium chloride
------------	---

関連情報

機能	May regulate calcium-dependent activities in the endoplasmic reticulum lumen or post-ER compartment.
配列類似性	Belongs to the CREC family. Contains 6 EF-hand domains.
細胞内局在	Endoplasmic reticulum lumen.

画像



15% SDS-PAGE analysis of 3µg ab98241

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors