

Recombinant Human RELM beta protein ab63292

画像数 2

製品の詳細

製品名	Recombinant Human RELM beta protein
精製度	> 95 % SDS-PAGE. Please filter the product by an appropriate sterile filter before using it in cell culture.
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MGSTQCSLDS VMDKKIKDVL NSLEYSPSPI SKKLSCASVK SQGRPSSCPA GMAVTGCACG YGCGSWDVQL ETTCHCQCSV VDWTTARCCH LT- KLRSHHHH HH

特性

Our **Abpromise guarantee** covers the use of **ab63292** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot SDS-PAGE
----------	--------------------------

製品の状態	Lyophilized
-------	-------------

前処理および保存


保存方法および安定性	Shipped at 4°C. Upon delivery aliquot. Store at -80°C. Avoid freeze / thaw cycle. Constituent: 0.41% Sodium acetate
------------	--

再構成	Add 0.1M Acetate buffer pH4 to prepare a working stock solution of approximately 0.5 mg/mL and let the lyophilized pellet dissolve completely. For conversion into higher pH value, we recommend intensive dilution by relevant buffer to a concentration of 10µg/mL. In higher concentrations the solubility of this antigen is limited.
-----	---


関連情報

機能	Probable hormone.
組織特異性	Expressed only in the gastrointestinal tract, particularly the colon.
配列類似性	Belongs to the resistin/FIZZ family.
細胞内局在	Secreted.

画像



LANE 1. M.W. marker – 14, 21, 31, 45, 66, 97 kDa. LANE 2. reduced and heated sample, 5µg/lane. LANE 3. non-reduced and non-heated sample, 5µg/lane.



Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors