

## Product datasheet

# Recombinant Human RHOD protein ab105608

画像数 1

### 製品の詳細

---

製品名	Recombinant Human RHOD protein
精製度	> 90 % SDS-PAGE. ab105608 is purified using conventional chromatography techniques.
発現系	Escherichia coli
アクセッション番号	<b><u>O00212</u></b>
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
配列	MMGSSHHHHHH SSGLVPRGSH MVKVVLVGDG GCGKTSLLMV FADGAFPESY TPTVFERYSMV NLQVKGKPVH LHIWDTAGQD DYDRLRPLFY PDASVLLLCF DVTSPNSFDN IFNRWYPEVN HFCKKVPIIV VGCKTDLRKD KSLVNKLRRN GLEPVTYHRG QEMARSVGAV AYLECSARLH DNVHAVFQEA AEVALSSRGR NFWRRITQGF C
予測される分子量	24 kDa including tags
領域	18 to 207
タグ	His tag N-Terminus

### 特性

---

Our **Abpromise guarantee** covers the use of **ab105608** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE Mass Spectrometry
質量分析	MALDI-TOF
製品の状態	Liquid

### 前処理および保存

---

## 保存方法および安定性

Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

pH: 8.00

Constituents: 0.00174% PMSF, 0.077% DTT, 0.316% Tris HCl, 0.0584% EDTA, 50% Glycerol (glycerin, glycerine), 1.16% Sodium chloride

## 関連情報

### 機能

Involved in endosome dynamics. May coordinate membrane transport with the function of the cytoskeleton. Participates in reorganization of actin cytoskeleton.

### 組織特異性

Heart, placenta, liver, skeletal muscle, and pancreas and, with weaker intensity, in several other tissues.


### 配列類似性

Belongs to the small GTPase superfamily. Rho family.

### 細胞内局在

Cell membrane.

## 画像



15% SDS-PAGE showing ab105608 (3 µg) at approximately 23.8 kDa.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors