

Recombinant Human RIDA protein ab101632

画像数 1

製品の詳細

製品名	Recombinant Human RIDA protein
精製度	> 95 % SDS-PAGE. ab101632 was purified using conventional chromatography.
発現系	Escherichia coli
アクセッション番号	<u>P52758</u>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	<p>MGSSHHHHHHSSGLVPRGSHMSSLIRRVISTAKAPGAIGP YSQAVLVDRT IYISGQIGMDPSSGQLVSGGVAEEAKQALKNMGEILKAAGCD FTNVVKTT VLLADINDFNTVNEIYKQYFKSNFPARAAYQVAALPKGSRIE IEAVAIQG PLTTASL</p>
予測される分子量	17 kDa including tags
領域	1 to 137
タグ	His tag N-Terminus

特性

Our **Abpromise guarantee** covers the use of **ab101632** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
	Mass Spectrometry
質量分析	MALDI-TOF
製品の状態	Liquid
備考	This product was previously labelled as HRSP12

前処理および保存

保存方法および安定性

Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

pH: 8.00

Constituents: 0.0154% DTT, 0.316% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

関連情報

機能

Endoribonuclease responsible for the inhibition of the translation by cleaving mRNA. Inhibits cell-free protein synthesis. Cleaves phosphodiester bonds only in single-stranded RNA.

組織特異性

Hepatocytes and renal distal tubular epithelial cells. Only weak expression in other tissues.

配列類似性

Belongs to the RutC family.


発生段階

Up-regulated during cellular differentiation.

細胞内局在

Cytoplasm. Nucleus. Mostly cytoplasmic but, in less differentiated cells occasionally nuclear.

画像



15% SDS PAGE analysis of 3 µg ab101632.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors