

## Recombinant human RIP3 protein (Tagged) ab125566

2 References [画像数 5](#)

### カルタヘナ法

#### 製品の詳細

製品名	Recombinant human RIP3 protein (Tagged)
生理活性	The specific activity of ab125566 was determined to be 38nmol/min/mg.
精製度	> 70 % Densitometry. Purity was determined to be >70% by densitometry.
発現系	Baculovirus infected Sf9 cells
アクセッション番号	<b><u>Q9Y572</u></b>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
予測される分子量	96 kDa including tags
領域	1 to 518
タグ	GST tag N-Terminus

#### 特性

Our **Abpromise guarantee** covers the use of **ab125566** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot Functional Studies SDS-PAGE
製品の状態	Liquid
備考	<b><u>ab64311</u></b> (Myelin Basic Protein protein) can be utilized as a substrate for assessing kinase activity

#### 前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 7.50 Constituents: 0.31% Glutathione, 0.002% PMSF, 0.004% DTT, 0.79% Tris HCl, 0.003% EDTA,
------------	---


25% Glycerol, 0.88% Sodium chloride

This product is an active protein and may elicit a biological response in vivo, handle with caution.


## 関連情報

機能	Promotes apoptosis.
組織特異性	Detected at lower levels in heart, placenta, lung and kidney. Isoform 3 is significantly increased in colon and lung cancers.
配列類似性	Belongs to the protein kinase superfamily. TKL Ser/Thr protein kinase family. Contains 1 protein kinase domain.
翻訳後修飾	Autophosphorylated.
細胞内局在	Cytoplasm.


## 画像



The specific activity of RIP3 (ab125566) was determined to be 10.6 nmol/min/mg as per activity assay protocol and was equivalent to 32 nmol/min/mg as per radiometric assay




SDS PAGE analysis of ab125566




SDS PAGE analysis of ab125566

SDS-PAGE - Recombinant human RIP3 protein  
(Tagged) (ab125566)



SDS-PAGE analysis of ab125566.

SDS-PAGE - Recombinant human RIP3 protein  
(Tagged) (ab125566)



The specific activity of ab125566 was determined to be 38 nmol/min/mg.

Functional Studies - Recombinant human RIP3  
protein (Tagged) (ab125566)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors