

## Product datasheet

# Recombinant Human RPS3 protein ab113158

画像数 1

### 製品の詳細

---

製品名	Recombinant Human RPS3 protein
精製度	> 90 % SDS-PAGE. ab113158 was purified by using conventional chromatography techniques.
発現系	Escherichia coli
アクセッション番号	<b><u>P23396</u></b>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	<b>MGSSHHHHHSSGLVPRGSHMAVQISKRRKFVADGIFKAE LNEFLTRELA EDGYSGVEVRVTPTRTEIIILATRTQNVLGEKRRIRELTAV VQKRFGFP EGSVELYAEKVATRGLCAIAQAESLRYKLLGGLAVRRACYGV LRFIMESG AKGCEVVVSGKLRGQRAKSMKFVDGLMIHSGDPVNYVDTAV RHVLLRQG VLGIKVKIMLPWDPTGKIGPKKPLPDHVSIVEPKDEILPTTP ISEQKGGK PEPPAMPQPVPTA</b>
予測される分子量	29 kDa
領域	1 to 263
タグ	His tag N-Terminus

### 特性

---

Our **Abpromise guarantee** covers the use of **ab113158** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE Mass Spectrometry
質量分析	MALDI-TOF
製品の状態	Liquid

## 前処理および保存

### 保存方法および安定性

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 0.02% DTT, 0.32% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

## 関連情報

### 配列類似性

Belongs to the ribosomal protein S3P family.

Contains 1 KH type-2 domain.

### 細胞内局在

Cytoplasm. Localized in cytoplasmic mRNP granules containing untranslated mRNAs.

## 画像



ab113158 at 3 µg analysed by 15% SDS PAGE. The molecular weight on SDS-PAGE will appear higher than predicted.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors