

Recombinant Human S100G protein ab96001

画像数 1

製品の詳細

製品名	Recombinant Human S100G protein
精製度	> 95 % SDS-PAGE. ab96001 is not sterile. Please filter the product by an appropriate sterile filter before using it in the cell culture.
発現系	Escherichia coli
アクセッション番号	<u>P29377</u>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MKHHHHHHASTKKSPPEELKRIFEKYAAKEGDPDQLSKDELK LLIQAEFPS LLKGPNTLDDLQELDKNGDGEVSFEFQVLVKKISQ
予測される分子量	10 kDa including tags
領域	2 to 79
タグ	His tag N-Terminus

特性

Our **Abpromise guarantee** covers the use of **ab96001** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot SDS-PAGE
製品の状態	Lyophilized

前処理および保存

保存方法および安定性	Shipped at 4°C. Upon delivery aliquot. Store at -80°C. Avoid freeze / thaw cycle. Constituents: 0.242% Tris, 0.29% Sodium chloride
再構成	Add deionized water to prepare a working stock solution of approximately 0.5 mg/mL and let the


lyophilized pellet dissolve completely.

関連情報

配列類似性

Belongs to the S-100 family.
Contains 2 EF-hand domains.

画像



14% SDS-PAGE separation of Human S100G protein (ab96001)

Lane 1. Molecular weight marker: 14, 21, 31, 45, 66 and 97 kDa

Lane 2. Reduced and heated sample, 5 µg per lane

Lane 3. Non reduced and non heated sample, 5 µg per lane

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors