

Recombinant Human S100P protein ab95488

1 References [画像数 1](#)

製品の詳細

製品名	Recombinant Human S100P protein
精製度	> 90 % SDS-PAGE. ab95488 was purified using conventional chromatography techniques
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MGSSHHHHHH SSGLVPRGSH MTELETAMGM IIDVFSRYSG SEGSTQTLTK GELKVLMEKE LPGFLQSGKD KDAVDKLLKD LDANGDAQVDFSEFIVFVAA ITSACHKYFE KAGLK
領域	1 to 95

特性

Our **Abpromise guarantee** covers the use of **ab95488** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE Mass Spectrometry ELISA
質量分析	MALDI-TOF-TOF
製品の状態	Liquid
備考	Endotoxin level is 1.0 EU per 1 microgram of protein (determined by LAL method).

前処理および保存


保存方法および安定性	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. pH: 8.00
------------	--

Constituents: 0.0154% DTT, 0.316% Tris HCl, 20% Glycerol (glycerin, glycerine), 0.29% Sodium chloride

関連情報

組織特異性	Up-regulated in various pancreatic ductal adenocarcinomas and pancreatic intraepithelial neoplasias.
配列類似性	Belongs to the S-100 family. Contains 2 EF-hand domains.
細胞内局在	Nucleus. Cytoplasm. Colocalizes with S100PBP in the nucleus.

画像



15% SDS-PAGE analysis of ab95488 (3ug).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors