

Recombinant Human SCGF alpha protein ab54145

製品の詳細

製品名	Recombinant Human SCGF alpha protein
精製度	> 95 % SDS-PAGE. Endotoxin level: < 0.1 ng per ug of SCGF alpha
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human

特性

Our **Abpromise guarantee** covers the use of **ab54145** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Lyophilized

前処理および保存

保存方法および安定性	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
再構成	The lyophilized SCGF alpha is soluble in water and most aqueous buffers. The lyophilized protein should be reconstituted in water to a concentration of 0.1-1.0 mg/ml. It is recommended that further dilutions be made into buffer containing carrier protein or medium containing serum. For most in vitro applications, SCGF alpha exerts its biological activity in the concentration range of 0.2 to 20.0 ng/ml.

関連情報

関連性	SCGF alpha (Stem Cell Growth Factor-alpha) is a newly discovered hematopoietic growth factor that exerts its activity at the early stages of hematopoiesis. Human SCGF alpha is a non-glycosylated 33.9 kDa species specific cytokine that can support growth of primitive hematopoietic cells, and in combination of EPO or GMCSF, it promotes proliferation of erythroid
-----	---

and myeloid progenitors, respectively.

細胞内局在

Cytoplasmic and Secreted

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors