

Recombinant Human SDSL protein ab113400

画像数 1

製品の詳細

製品名	Recombinant Human SDSL protein
精製度	> 90 % SDS-PAGE. Conventional chromatography techniques (anion exchange followed by gel filtration with 20mM Tris pH 7.5, 2mM EDTA).
発現系	Escherichia coli
アクセッション番号	<u>Q96GA7</u>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	<p>MGSSHHHHHH SSGLVPRGSH MGSHMDGPVA EHAKQEPFHV VTPLLESWAL SQVAGMPVFL KCENVQPSGS FKIRGIGHFC QEMAKKGRHLVCSSGGNAG IAAAYAARKL GIPATIVLPE STSLQVVQRL QGEGAEVQLT GKVDDEANLR AQELAKRDGW ENVPPFDHPL IWKGHASLVQ MGSHMDGPVA EHAKQEPFHV VTPLLESWAL SQVAGMPVFL KCENVQPSGS FKIRGIGHFC QEMAKKGRH ELKAVLRTPP GALVLAVGGG GLLAGVVAGL LEVGWQHVPI IAMETHGAHC FNAAITAGKL VTLDPDITSVA KSLGAKTVAA RALECMQVCK IHSEVVEDTE AVSAVQQLLD DERMLVEPAC GAALAAIYSG LLRRQLAEGC LPPSLTSVVV IVCGGNNINS RELQALKTHL GQV</p>
予測される分子量	37 kDa including tags
領域	1 to 329
タグ	His tag N-Terminus

特性

Our **Abpromise guarantee** covers the use of **ab113400** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Mass Spectrometry
	SDS-PAGE
質量分析	MALDI-TOF
製品の状態	Liquid


前処理および保存

保存方法および安定性	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. pH: 8.00 Constituents: 0.04% DTT, 0.32% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.58% Sodium chloride
------------	--

関連情報

機能	Has low serine dehydratase and threonine dehydratase activity.
配列類似性	Belongs to the serine/threonine dehydratase family.

画像



15% SDS-PAGE of SDSL protein: (1) MW ladder, (2) protein loaded 3ug

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors