

Recombinant Human SF20 protein ab86915

1 References [画像数 1](#)

製品の詳細

製品名	Recombinant Human SF20 protein
精製度	> 95 % SDS-PAGE. ab86915 was purified using conventional chromatography techniques, after refolding of the isolated inclusion bodies in a renaturation buffer.
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	<p>MGSSHHHHHH SSGLVPRGSH MSEPTTVAFD</p> <p>VRPGGVVHSF SHNVGPGDKY TCMFTYASQG</p> <p>GTNEQWQMSL GTSEDHQHFT CTIWRPQGKS</p> <p>YLYFTQFKAE VRGAEIEYAM AYSKAAFERE</p> <p>SDVPLKTEEF EVTKTAVAGR PGAFKAELSK</p> <p>LVIVAKASRT EL</p>

特性

Our **Abpromise guarantee** covers the use of **ab86915** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション SDS-PAGE

製品の状態 Liquid

前処理および保存

保存方法および安定性 Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

pH: 8.00

Constituents: 0.316% Tris HCl, 20% Glycerol (glycerin, glycerine)

関連情報

組織特異性	Expressed in synovial tissue. Found in synovial fluid of patients with arthropaties.
配列類似性	Belongs to the UPF0556 family.
細胞内局在	Secreted. Endoplasmic reticulum-Golgi intermediate compartment.

画像



15% SDS-PAGE analysis of 3 ug ab86915

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors