

### Recombinant Human Slit3 protein ab116753

画像数 1

#### 製品の詳細

製品名	Recombinant Human Slit3 protein
発現系	Wheat germ
アクセッション番号	<b><u>O75094</u></b>
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
配列	PCLGHRCHHGKCVATGTSYMCKCAEGYGGDLCDNKNSANAC SAFKCHHG QCHISDQGEPTYCLCQPGFSGEHCQQENPCLGQVVREVIRRQK GYASCATA
予測される分子量	37 kDa including tags
領域	1371 to 1470

#### 特性

Our **Abpromise guarantee** covers the use of **ab116753** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	ELISA
	SDS-PAGE
	Western blot
製品の状態	Liquid


#### 前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.3% Glutathione, 0.79% Tris HCl
------------	---

#### 関連情報

機能	May act as molecular guidance cue in cellular migration, and function may be mediated by interaction with roundabout homolog receptors.
組織特異性	Predominantly expressed in thyroid.
配列類似性	<p>Contains 1 CTCK (C-terminal cystine knot-like) domain.</p> <p>Contains 9 EGF-like domains.</p> <p>Contains 1 laminin G-like domain.</p> <p>Contains 20 LRR (leucine-rich) repeats.</p> <p>Contains 4 LRRCT domains.</p> <p>Contains 4 LRRNT domains.</p>
細胞内局在	Secreted.

## 画像



12.5% SDS-PAGE analysis of ab116753 at approximately 36.63kDa, stained with Coomassie Blue.

SDS-PAGE - Recombinant Human Slit3 protein  
(ab116753)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors