

Recombinant Human SMARCA6 protein ab114875

画像数 1

製品の詳細

製品名	Recombinant Human SMARCA6 protein
発現系	Wheat germ
アクセッション番号	<u>Q9NRZ9</u>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	<p>MFGSSEKETIELSPTGRPKRRTRKSINYSKIDDFPNELEKLI SQIQPEVD RERAVVEVNIPVESEVNLKLQNIMMLLRKCCNHPYLIEYPID PVTQEFKI DEELVTNSGKFLILDRMLPELKKRGHKVLLFSQMTSMLDILM DYCHLRDF NFSRLDGSMYSEREKNMHSFNTDPEVFIFLVSTRAGGLGIN LTAADTVI IYDSDWNPQSDLQAQDRCHRIGQTKPVVVYRLVTANTIDQKI VERAAAKR KLEKLIHKNHFKGGQSGNLNLSKNFLDPKELMELLKSRDYER EIKGSREK VISDKDLELLLDLRSIDLIDQMNASGPIKEKMGIFKILENSED SPECLF</p>
予測される分子量	64 kDa including tags
領域	1 to 348

特性

Our **Abpromise guarantee** covers the use of **ab114875** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	ELISA
	SDS-PAGE
	Western blot
製品の状態	Liquid

前処理および保存

保存方法および安定性

Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.3% Glutathione, 0.79% Tris HCl

関連情報

機能

Plays an essential role in normal development and survival. Involved in regulation of the expansion or survival of lymphoid cells. Required for de novo or maintenance DNA methylation. May control silencing of the imprinted CDKN1C gene through DNA methylation. May play a role in formation and organization of heterochromatin, implying a functional role in the regulation of transcription and mitosis.

組織特異性

Highly expressed in proliferative tissues such as adult thymus and testis, and expressed at lower levels in uterus, small intestine, colon, and peripheral blood mononuclear cells. Also expressed in neoplastic cell lines including those derived from myeloid and lymphoid leukemias.

配列類似性

Belongs to the SNF2/RAD54 helicase family.


Contains 1 helicase ATP-binding domain.

Contains 1 helicase C-terminal domain.

細胞内局在

Nucleus. Closely associated with pericentric heterochromatin.

画像



12.5% SDS-PAGE showing ab114875 at approximately 64.39kDa.

Stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors