

### Recombinant Human SOX2 protein ab79950

画像数 2

#### 製品の詳細

製品名	Recombinant Human SOX2 protein
精製度	> 95 % SDS-PAGE. ab79950 is greater than 95% by SDS-PAGE gel and HPLC analyses. Endotoxin level is less than 0.1 ng per µg (1EU/µg).
エンドキシン・レベル	< 0.100 Eu/µg
発現系	Escherichia coli
アクセッション番号	<b><u>P48431</u></b>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human

配列	MYNMMETELK PPGPQQTSGG GGGNSTAAAA GGNQKNPDR VKRPMNAFMV WSRGQRRKMA QENPKMHNSE ISKRLGAEWK LLSETEKRPF IDEAKRLRAL HMKEHPDYKY RPRRKTCTLM KKDKYTLPGG LLAPGGNSMA SGVGVGAGLG AGVNQRMDSY AHMNGWSNGS YSMMQDQLGY PQHPGLNAHG AAQMMPMHRY DVSALQYNM TSSQTYMNGS PTYSMSYSQQ GTPGMALGSM GSVVKSEASS SPPVVTSSSH SRAPCQAGDL RDMISMYLPG AEVPEPAAPS RLHMSQHYQS GPVPGTAING TLPLSHM
----	--

領域	1 to 317
----	----------

#### 特性

Our **Abpromise guarantee** covers the use of **ab79950** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Sandwich ELISA
	Western blot
	HPLC
	SDS-PAGE

## 製品の状態

Lyophilized

## 前処理および保存

### 保存方法および安定性

Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.

### 再構成

Reconstituted ab79950 is stable for at least 3 months when stored in working aliquots with a carrier protein at -20°C. Avoid repeated freeze/thaw cycles.

## 関連情報

### 機能

Transcription factor that forms a trimeric complex with OCT4 on DNA and controls the expression of a number of genes involved in embryonic development such as YES1, FGF4, UTF1 and ZFP206 (By similarity). Critical for early embryogenesis and for embryonic stem cell pluripotency.

### 関連疾患

Defects in SOX2 are the cause of microphthalmia syndromic type 3 (MCOPS3) [MIM:206900]. Microphthalmia is a clinically heterogeneous disorder of eye formation, ranging from small size of a single eye to complete bilateral absence of ocular tissues (anophthalmia). In many cases, microphthalmia/anophthalmia occurs in association with syndromes that include non-ocular abnormalities. MCOPS3 is characterized by the rare association of malformations including uni- or bilateral anophthalmia or microphthalmia, and esophageal atresia with trachoesophageal fistula.

### 配列類似性

Contains 1 HMG box DNA-binding domain.


### 翻訳後修飾

Sumoylation inhibits binding on DNA and negatively regulates the FGF4 transactivation.


### 細胞内局在

Nucleus.

## 画像



Standard Curve for SOX2 (Analyte: **SOX2 protein (Human) (ab79950)**); dilution range 1pg/ml to 1µg/ml using Capture Antibody **Mouse monoclonal [57CT23.3.4] to SOX2 (ab75485)** at 0.2µg/ml and Detector Antibody **Rabbit monoclonal [EPR3131] to SOX2 (ab92494)** at 0.5µg/ml.



**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors