

# Recombinant Human Stanniocalcin 1/STC protein ab63280

画像数 1

### 製品の詳細

製品名	Recombinant Human Stanniocalcin 1/STC protein
精製度	> 95 % SDS-PAGE. Please filter the product by an appropriate sterile filter before using it in the cell culture.
発現系	HEK 293 cells
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	THEAEQNDSV SPRKSRVAAQ NSAEVVRCLN SALQVGC GAF ACLENSTCDT DGMVDICKSF LYSAAKFDTQ GKAFVKESLK CIANGVTSKV FLAIRRCSTF QRMIAEVQEE CYSKLVCSI AKRNPEAITE VVQLPNHFSN RYYNRLVRS LECEDTVST IRDSLMEKIG PNMASLFHIL QTDHCAQTHP RADFNRRRTN EPQKLKVLRL NLRGEEDSPS HIKRTSHESA - ASDYKDDDDK

### 特性

Our **Abpromise guarantee** covers the use of **ab63280** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot ELISA SDS-PAGE
製品の状態	Lyophilized
備考	Previously labelled as Stanniocalcin 1.

### 前処理および保存

保存方法および安定性	Shipped at 4°C. Upon delivery aliquot. Store at -80°C. Avoid freeze / thaw cycle. Constituents: 0.242% Tris, 0.116% Sodium chloride
------------	--

**再構成** Add deionized water to prepare a working stock solution of approximately 0.5 mg/mL and let the lyophilized pellet dissolve completely.

## 関連情報


**機能** Stimulates renal phosphate reabsorption, and could therefore prevent hypercalcemia.

**組織特異性** Expressed in most tissues, with the highest levels in ovary, prostate, heart, kidney and thyroid. In the kidney, expression is confined to the nephron, specifically in the distal convoluted tubule and in the collecting tubule. Not detected in the brain, liver, spleen, peripheral blood leukocytes and adrenal medulla.

**配列類似性** Belongs to the stanniocalcin family.

**細胞内局在** Secreted.

## 画像



LANE 1. M.W. marker - 14, 21, 31, 45, 66, 97 kDa. LANE 2. reduced and heated sample, 5µg/lane. LANE 3. non-reduced and non-heated sample, 5µg/lane

SDS-PAGE - Recombinant Human Stanniocalcin 1/STC protein (ab63280)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors