

# Recombinant Human STIP1/STI1 protein ab92357

画像数 2

### 製品の詳細

---

製品名	Recombinant Human STIP1/STI1 protein
精製度	> 90 % SDS-PAGE. >90% pure as determined by SDS-PAGE and Western blot analyses. Purification by multi step chromatography.
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human

### 特性

---

Our **Abpromise guarantee** covers the use of **ab92357** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE Western blot
製品の状態	Liquid
備考	This product was previously labelled as STIP1

### 前処理および保存

---

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. Constituent: PBS
------------	---

### 関連情報

---

機能	Mediates the association of the molecular chaperones HSC70 and HSP90 (HSPCA and HSPCB).
配列類似性	Contains 2 STI1 domains.

Contains 9 TPR repeats.


## ドメイン

The TPR 1 repeat interacts with the C-terminal of HSC70. The TPR 4, 5 and 6 repeats (also called TPR2A domain) and TPR 7, 8 and 9 repeats (also called TPR2B domain) interact with HSP90.

## 細胞内局在

Cytoplasm. Nucleus.

## 画像



SDS-PAGE - Recombinant Human STIP1/STI1 protein (ab92357)

SDS-PAGE analysis of ab92357 (STIP1/STI1 protein) stained with Coomassie Blue.


Lane 1: Molecular Weight marker

Lane 2: 0.5 µg

Lane 3: 1 µg

Lane 4: 2 µg

Lane 5: 5 µg



Western blot - Recombinant Human STIP1/STI1 protein (ab92357)

**All lanes :** STIP1/STI1 antibody at 1 µg/ml

**Lane 1 :** Molecular weight marker

**Lane 2 :** STIP1/STI1 protein (ab92357) at 0.05 µg

**Lane 3 :** STIP1/STI1 protein (ab92357) at 0.1 µg

**Predicted band size:** 63 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors