

Recombinant Human STK38 protein ab101505

1 References [画像数 1](#)

製品の詳細

製品名	Recombinant Human STK38 protein
精製度	> 80 % Densitometry. The purity was determined to be >80% by densitometry. Affinity purified.
発現系	Baculovirus infected insect cells
アクセッション番号	2.7.11.1
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
予測される分子量	81 kDa including tags
領域	2 to 465

特性

Our **Abpromise guarantee** covers the use of **ab101505** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot SDS-PAGE
製品の状態	Liquid

前処理および保存


保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 7.50 Constituents: 0.307% Glutathione, 0.00174% PMSF, 0.00385% DTT, 0.79% Tris HCl, 0.00292% EDTA, 25% Glycerol (glycerin, glycerine), 0.87% Sodium chloride
------------	--

関連情報

機能	Negative regulator of MAP3K1/2 signaling. Converts MAP3K2 from its phosphorylated form to its
----	---

	nonphosphorylated form and inhibits autophosphorylation of MAP3K2.
組織特異性	Ubiquitously expressed with highest levels observed in peripheral blood leukocytes.
配列類似性	Belongs to the protein kinase superfamily. AGC Ser/Thr protein kinase family. Contains 1 AGC-kinase C-terminal domain. Contains 1 protein kinase domain.
翻訳後修飾	ISGylated.
細胞内局在	Nucleus. Cytoplasm.

画像



SDS-PAGE analysis of ab101505.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors