

### Recombinant Human STT3A protein ab116946

画像数 1

#### 製品の詳細

製品名	Recombinant Human STT3A protein
発現系	Wheat germ
アクセッション番号	<b><u>P46977</u></b>
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
配列	GGSTDTGKHIKENDYYTPTGFEFRVDREGSPVLLNCLMYKMCY YRFGQVYT EAKRPPGFDRVRNAEIGNKDFELDVLEEAYTTEHWLVRIYKV KDLDNRG
予測される分子量	37 kDa including tags
領域	603 to 701

#### 特性

Our **Abpromise guarantee** covers the use of **ab116946** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	ELISA
	Western blot
	SDS-PAGE
製品の状態	Liquid


#### 前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.3% Glutathione, 0.79% Tris HCl
------------	---

#### 関連情報

<b>機能</b>	Component of the N-oligosaccharyl transferase enzyme which catalyzes the transfer of a high mannose oligosaccharide from a lipid-linked oligosaccharide donor to an asparagine residue within an Asn-X-Ser/Thr consensus motif in nascent polypeptide chains. N-glycosylation occurs cotranslationally and the complex associates with the Sec61 complex at the channel-forming translocon complex that mediates protein translocation across the endoplasmic reticulum (ER). SST3A seems to be involved in complex substrate specificity.
<b>組織特異性</b>	Expressed at high levels in placenta, liver, muscle and pancreas, and at very low levels in brain, lung and kidney. Expressed in skin fibroblasts (at protein level).
<b>配列類似性</b>	Belongs to the STT3 family.
<b>細胞内局在</b>	Endoplasmic reticulum membrane.

## 画像



12.5% SDS-PAGE showing ab116946 at approximately 36.52kDa stained with Coomassie Blue.

SDS-PAGE - Recombinant Human STT3A protein  
(ab116946)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors