

Recombinant Human TCEAL1 protein ab108377

画像数 1

製品の詳細

製品名	Recombinant Human TCEAL1 protein
精製度	> 90 % SDS-PAGE. ab108377 was purified using conventional chromatography techniques.
発現系	Escherichia coli
アクセッション番号	<u>Q15170-2</u>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MDKPRKENEEEPQSAPKTDEERPPVEHSPEKQSPEEQSSEEQ SSEEEFFP EELLPELLPEMLLSEERPPQEGLSRKDLFEGRPPMEQPPCGV GKHKLEEG SFKERLARSRPQFRGDIHGRNLSNEEMIQAADLEEMKVRN KLMIMHWK AKRSRPYPPILEHHHHHH
予測される分子量	20 kDa including tags
領域	1 to 159
タグ	His tag C-Terminus

特性

Our **Abpromise guarantee** covers the use of **ab108377** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE Mass Spectrometry
質量分析	MALDI-TOF
製品の状態	Liquid

前処理および保存

保存方法および安定性	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -
------------	--

80°C. Avoid freeze / thaw cycle.


pH: 8.00

Constituents: 0.0308% DTT, 0.316% Tris HCl, 20% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

関連情報

機能	May be involved in transcriptional regulation. Modulates various viral and cellular promoters in a promoter context-dependent manner. For example, transcription from the FOS promoter is increased, while Rous sarcoma virus (RSV) long terminal repeat (LTR) promoter activity is repressed. Does not bind DNA directly.
組織特異性	Expressed in all tissues examined. Highly expressed in heart, ovary, prostate and skeletal muscle. Moderately expressed in brain, placenta, testis and small intestine. Weakly expressed in lung, liver and spleen. Expressed in several cancer cell lines.
配列類似性	Belongs to the TFS-II family. TFA subfamily.
翻訳後修飾	Phosphorylation of Ser-36 and Ser-37 is critical for transcriptional repression.
細胞内局在	Nucleus.

画像



15% SDS-PAGE analysis of ab108377 (3µg)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors