

Recombinant Human TNNI1 protein ab125609

画像数 1

製品の詳細

製品名	Recombinant Human TNNI1 protein
精製度	> 75 % Densitometry. Purity determined to be >75% by densitometry.
発現系	Escherichia coli
アクセッション番号	<u>P19237</u>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
予測される分子量	24 kDa including tags
領域	1 to 187
タグ	His tag N-Terminus

特性

Our **Abpromise guarantee** covers the use of **ab125609** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot SDS-PAGE
製品の状態	Liquid


前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 7.00 Preservative: 1.02% Imidazole Constituents: 0.002% PMSF, 0.81% Sodium phosphate, 0.03% DTT, 25% Glycerol (glycerin, glycerine), 1.75% Sodium chloride
------------	---

関連情報

機能	Troponin I is the inhibitory subunit of troponin, the thin filament regulatory complex which confers calcium-sensitivity to striated muscle actomyosin ATPase activity.
組織特異性	Highest levels observed in human skeletal muscle (e.g. gastrocnemius muscle), differentiated cultures of primary human muscle cells and rhabdomyosarcoma cells cultured in low serum medium. Expressed in C2 muscle cell myoblasts and myotubes.
配列類似性	Belongs to the troponin I family.

画像



SDS-PAGE analysis of ab125609.

SDS-PAGE - Recombinant Human TNNI1 protein
(ab125609)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors