

Product datasheet

Recombinant Human TRA2B/SFRS10 protein ab114393

画像数 1

製品の詳細

製品名	Recombinant Human TRA2B/SFRS10 protein
発現系	Escherichia coli
アクセッション番号	P62995
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
配列	LGVFGLSLYTTERDLREVFSKYGPIADVSIYDQQSRRSRGF AFVYFENV DDAKEAKERANGMELDGRRIRVDFSITKRP
予測される分子量	34 kDa including tags
領域	120 to 199

特性

Our **Abpromise guarantee** covers the use of **ab114393** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE Western blot ELISA
製品の状態	Liquid
備考	This product was previously labelled as TRA2B.


前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.79% Tris citrate/phosphate, 0.3% Glutathione
------------	---

関連情報

機能	Sequence-specific RNA-binding protein which participates in the control of pre-mRNA splicing.
組織特異性	Highest expression in heart, skeletal muscle and pancreas. Less abundant in kidney, placenta and brain. Lowest expression in kidney and liver.
配列類似性	Belongs to the splicing factor SR family. Contains 1 RRM (RNA recognition motif) domain.
翻訳後修飾	Phosphorylated in the RS domains. Dimethylation at Arg-241 is probably asymmetric.
細胞内局在	Nucleus.

画像



ab114393 analysed on a 12.5% SDS-PAGE Stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors