

# Recombinant Human Urokinase protein ab92591

画像数 1

### 製品の詳細

製品名	Recombinant Human Urokinase protein
精製度	> 95 % SDS-PAGE. ab92591 is purified from auto-catalytic digestion of recombinant human urokinase produced in insect cells.
発現系	Insect cells
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
領域	1 to 135

### 特性

Our **Abpromise guarantee** covers the use of **ab92591** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Liquid
備考	The amino terminal fragment (ATF) contains the Kringle and Epidermal Growth Factor domains and is separated from LMW uPA by cleavage at the Lys135-Lys136 bond. Produced using non-baculovirus insect cells.


### 前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 7.40 Constituents: 2.38% HEPES, 0.58% Sodium chloride
------------	---

### 関連情報

機能	Specifically cleave the zymogen plasminogen to form the active enzyme plasmin.
組織特異性	Expressed in the prostate gland and prostate cancers.
配列類似性	Belongs to the peptidase S1 family. Contains 1 EGF-like domain. Contains 1 kringle domain. Contains 1 peptidase S1 domain.
翻訳後修飾	Phosphorylation of Ser-158 and Ser-323 abolishes proadhesive ability but does not interfere with receptor binding.
細胞内局在	Secreted.

## 画像



### 12% SDS-PAGE

Lane 1: Urokinase (2 ug) Non-reduced

Lane 2: Prestained Standard

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors