

Product datasheet

Recombinant Human VAMP2 protein ab49001

画像数 1

製品の詳細

| | |
|-------------|---|
| 製品名 | Recombinant Human VAMP2 protein |
| 精製度 | > 95 % SDS-PAGE. |
| 発現系 | Escherichia coli |
| アクセッション番号 | <u>P63027</u> |
| タンパク質長 | Protein fragment |
| Animal free | No |
| 由来 | Recombinant |
| 生物種 | Human |
| 配列 | MRGSHHHHHH GMASMTGGQQ MGRDLYDDDD KDRWGSMSA TAATAPPAAP AGEGGPPAPP PNLTSNRRLQ QTQAQVDEVV DIMRVNVDKV LERDQKLESEL DDRADALQAG ASQFETSAAK LKRKYW |
| 予測される分子量 | 14 kDa |
| 領域 | 1 to 89 |

特性

Our **Abpromise guarantee** covers the use of **ab49001** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| | |
|----------|----------|
| アプリケーション | SDS-PAGE |
| 製品の状態 | Liquid |


前処理および保存

| | |
|------------|--|
| 保存方法および安定性 | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. pH: 7.40 Constituents: PBS, 0.0292% EDTA |
|------------|--|

関連情報

| | |
|-------|---|
| 機能 | Involved in the targeting and/or fusion of transport vesicles to their target membrane. |
| 組織特異性 | Nervous system and skeletal muscle. |
| 配列類似性 | Belongs to the synaptobrevin family. Contains 1 v-SNARE coiled-coil homology domain. |
| 細胞内局在 | Cytoplasmic vesicle > secretory vesicle > synaptic vesicle membrane. Cell junction > synapse > synaptosome. Neuronal synaptic vesicles. |

画像



ab49001 run on a 15% SDS-PAGE gel with molecular weight markers.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors