

### Recombinant Human Vaspin protein ab50136

#### 製品の詳細

製品名	Recombinant Human Vaspin protein
精製度	> 95 % SDS-PAGE. Greater than 98% by SDS-PAGE and HPLC analyses. Endotoxin level is less than 0.1 ng per µg (1EU/µg).
エンドトキシン・レベル	< 0.100 Eu/µg
発現系	Escherichia coli
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
配列	<div> MLKPSFSPRN YKALSEVQGW KQRMAAKELA  RQNMDLGFKL LKKLAFYNPG RNIFLSPLSI  STAFSMLCLG AQDSTLDEIK QGFNFRKMPE  KDLHEGFHYI IHELTQKTQD LKLSIGNTLF  IDQRLQPQRK FLEDAKNFYS AETILTNFQN  LEMAQKQIND FISQKTHGKI NNLIENIDPG  TVMLLANYIF FRARWKHEFD PNVTKEDFF  LEKNSSVKVP MMFRSGIYQV GYDDKLSCTI  LEIPYQKNIT AIFILPDEGK LKHLEKGLQV  DTFSRWKTLL SRRVVDVSVR RLHMTGTFDL  KKTLSYIGVS KIFEEHGDLT KIAPHRSLKV  GEAVHKAELK MDERGTEGAA GTGAQTL PME  TPLVVKIDKP YLLLIYSEKI PSVLFLGKIV NPIGK </div>

#### 特性

Our **Abpromise guarantee** covers the use of **ab50136** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Lyophilized
備考	This product was previously labelled as Vaspin, SERPINA12

## 前処理および保存

---

### 保存方法および安定性

Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

Constituents: 0.87% Sodium chloride, 0.61% Tris

### 再構成

Centrifuge the vial prior to opening. Reconstitute in water to a concentration of 0.1-1.0 mg/ml. This solution can then be diluted into other aqueous buffers and stored at 4°C for 1 week or -20°C for future use.

## 関連情報

---

機能	May modulates insulin action conceivably only in the presence of its yet undefined target proteases in white adipose tissues.
組織特異性	Expressed in visceral adipose tissues.
配列類似性	Belongs to the serpin family.
細胞内局在	Secreted.

---

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors