

Recombinant Human Vaspin protein ab86032

1 References [画像数 1](#)

製品の詳細

製品名	Recombinant Human Vaspin protein
精製度	> 90 % SDS-PAGE. ab86032 is purified using conventional chromatography techniques.
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MGSSHHHHHH SSGLVPRGSH MLKPSFSPRN YKALSEVQGW KQRMAAKELA RQNMDLGFKL LKKLAFYNPG RNIFLSPLSI STAFSMLCLG AQDSTLDEIK QGFNFRKMPE KDLHEGFHYI IHELTQKTQD LKLSIGNTLF IDQRLQPQRK FLEDAKNFYS AETILTNTFNQ LEMAKQKIND FISQKTHGKI NNLIENIDPG TVMLLANYIF FRARWKHEFD PNVTKKEEDFF LEKNSSVKVP MMFRSGIYQV GYDDKL SCTI LEIPYQKNIT AIFILPDEGK LKHLEKGLQV DTFSRWKTLL SRRVVDVSV R LHMTGTFDL KKTLSYIGVS KIFEEHGDLT KIAPHRSLKV GEAVHKAELK MDERGTEGAA GTGAQTL PME TPLVVKIDKP YLLLIYSEKI PSVLFLGKIV NPIGK

特性

Our **Abpromise guarantee** covers the use of **ab86032** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Liquid
備考	This product was previously labelled as Vaspin, SERPINA12

前処理および保存

保存方法および安定性

Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.


pH: 8.00

Constituents: 0.00348% PMSF, 0.316% Tris HCl, 10% Glycerol (glycerin, glycerine)

関連情報

機能	May modulates insulin action conceivably only in the presence of its yet undefined target proteases in white adipose tissues.
組織特異性	Expressed in visceral adipose tissues.
配列類似性	Belongs to the serpin family.
細胞内局在	Secreted.

画像



SDS-PAGE analysis of 3 µg ab86032 under reducing conditions and visualized by coomassie blue stain.

SDS-PAGE - Recombinant Human Vaspin protein
(ab86032)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors