

Product datasheet

Recombinant Human VP16 protein ab81832

★★★★★ 1 Abreviews

製品の詳細

製品名	Recombinant Human VP16 protein
精製度	> 95 % SDS-PAGE. Isolated from an E. coli strain that carries the coding sequence of the fused protein under the control of a T7 promoter.
発現系	Escherichia coli
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
配列	STAPPTDVSLGDELHLDGEDVAMAHADALDDFDLMLGDGDS PGPGFTPH DSAPYGALDMADFEFEQMFTDALGIDEYGG
領域	411 to 490
タグ	GST tag N-Terminus

特性

Our **Abpromise guarantee** covers the use of **ab81832** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Liquid

前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 7.9 Constituents: 0.75% Potassium chloride, 0.0154% DTT, 0.316% Tris HCl, 0.00584% EDTA, 20% Glycerol (glycerin, glycerine)
------------	---

関連情報

関連性

Herpes virus VP16 activates expression of immediate early genes in virally-infected cells. As most other eukaryotic transcriptional activator proteins, VP16 has a modular domain structure: it's N-terminus is involved in DNA-protein interactions, while it's C-terminal 79 amino acids have proven to be an especially potent transactivation domain. VP16 has been shown to bind to TBP, TFIIIB, and replication factor A

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors