

Recombinant human VRK2 protein ab107698

画像数 5

製品の詳細

製品名	Recombinant human VRK2 protein
生理活性	Specific activity: 12 nmol/min/mg
精製度	> 90 % SDS-PAGE. The purity was determined to be >90% by densitometry. Affinity purified.
発現系	Baculovirus infected Sf9 cells
アクセッション番号	<u>Q86Y07</u>
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
予測される分子量	66 kDa including tags
領域	1 to 375

特性

Our **Abpromise guarantee** covers the use of **ab107698** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE Functional Studies
製品の状態	Liquid
備考	<u>ab91092</u> (Cow Casein full length protein) can be utilized as a substrate for assessing Kinase activity


前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 7.50 Constituents: 0.307% Glutathione, 0.00174% PMSF, 0.00385% DTT, 0.79% Tris HCl, 0.00292% EDTA, 25% Glycerol (glycerin, glycerine), 0.87% Sodium chloride This product is an active protein and may elicit a biological response in vivo, handle with caution.
------------	--


関連情報

機能	Probable serine/threonine kinase.
組織特異性	Widely expressed. Highly expressed in fetal liver, skeletal muscle, pancreas, heart, peripheral blood leukocytes and testis.
配列類似性	Belongs to the protein kinase superfamily. CK1 Ser/Thr protein kinase family. VRK subfamily. Contains 1 protein kinase domain.
細胞内局在	Membrane.


画像



The specific activity of VRK2 (ab107698) was determined to be 10 nmol/min/mg as per activity assay protocol




SDS PAGE analysis of ab107698



SDS PAGE analysis of ab107698



SDS-PAGE analysis of ab107698



Specific activity of ab107698 was determined to be 12 nmol/min/mg

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors