

### Recombinant KIR2DL1 protein ab51212

画像数 1

#### 製品の詳細

製品名	Recombinant KIR2DL1 protein
精製度	> 95 % SDS-PAGE. ab51212 was purified by FPLC gel-filtration chromatography, after refolding of the isolated inclusion bodies in a redox buffer.
発現系	Escherichia coli
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant

#### 特性

Our **Abpromise guarantee** covers the use of **ab51212** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Liquid
備考	ab51212 corresponds to the extracellular domain of KIR2DL1 and was cloned into an E. coli expression vector and overexpressed as insoluble protein aggregates (inclusion bodies).


#### 前処理および保存

保存方法および安定性	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. pH: 7.50 Constituent: 0.316% Tris HCl
------------	---

#### 関連情報

機能	Receptor on natural killer (NK) cells for HLA-C alleles. Inhibits the activity of NK cells thus preventing cell lysis.
配列類似性	Belongs to the immunoglobulin superfamily. Contains 2 Ig-like C2-type (immunoglobulin-like) domains.

## 画像



14% SDS-PAGE gel loaded with ab51212, recombinant human KIR2DL1 protein (extracellular domain).

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

**Terms and conditions**

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors