

Recombinant KIR2DL3 protein ab51215

画像数 1

製品の詳細

| | |
|-------------|---|
| 製品名 | Recombinant KIR2DL3 protein |
| 精製度 | > 95 % SDS-PAGE. ab51215 was purified by FPLC gel-filtration chromatography. |
| 発現系 | Escherichia coli |
| タンパク質長 | Protein fragment |
| Animal free | No |
| 由来 | Recombinant |

特性

Our **Abpromise guarantee** covers the use of **ab51215** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| | |
|----------|----------|
| アプリケーション | SDS-PAGE |
| 製品の状態 | Liquid |


前処理および保存

| | |
|------------|---|
| 保存方法および安定性 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. pH: 7.50 Constituent: 0.316% Tris HCl |
|------------|---|

関連情報

| | |
|-------|--|
| 機能 | Receptor on natural killer (NK) cells for HLA-C alleles (HLA-Cw1, HLA-Cw3 and HLA-Cw7). Inhibits the activity of NK cells thus preventing cell lysis. |
| 配列類似性 | Belongs to the immunoglobulin superfamily. Contains 2 Ig-like C2-type (immunoglobulin-like) domains. |
| 細胞内局在 | Cell membrane. |

画像



14% SDS-PAGE gel loaded with ab51215, recombinant human KIR2DL3 protein (extracellular domain).
Predicted molecular weight 22.2 kDa

SDS-PAGE - Recombinant KIR2DL3 protein
(ab51215)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors