

Recombinant Mouse CXCL11 protein ab50170

製品の詳細

製品名	Recombinant Mouse CXCL11 protein
精製度	> 95 % SDS-PAGE. ab50170 purity was assessed also by HPLC. Endotoxin level is less than 0.1 ng per µg (1EU/µg).
エンドトキシン・レベル	< 0.100 Eu/µg
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Mouse
配列	FLMFKQGRCL CIGPGMKAVK MAEIEKASVI YPSNGCDKVE VIVTMKAHQR QRCLDPRSKQ ARLIMQAIEK KNFLRRQNM
領域	22 to 100

特性

Our **Abpromise guarantee** covers the use of **ab50170** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE Functional Studies
製品の状態	Lyophilized

前処理および保存

保存方法および安定性	Shipped at 4°C. The lyophilized protein is stable for a few weeks at room temperature. Store at -20°C long term.
再構成	Reconstitute to 1mg/ml in distilled water.

関連情報

機能	Chemotactic for interleukin-activated T-cells but not unstimulated T-cells, neutrophils or
----	--

	monocytes. Induces calcium release in activated T-cells. Binds to CXCR3. May play an important role in CNS diseases which involve T-cell recruitment. May play a role in skin immune responses.
組織特異性	High levels in peripheral blood leukocytes, pancreas and liver astrocytes. Moderate levels in thymus, spleen and lung. Low levels in placenta, prostate and small intestine. Also found in epidermal basal layer keratinocytes in skin disorders.
配列類似性	Belongs to the intercrine alpha (chemokine CxC) family.
細胞内局在	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors