

# Recombinant Mouse FGF21 + Human IgG1 fusion protein (Fc Chimera) ab108580

### 製品の詳細

|             |   |
|-------------|---|
| 製品名         | Recombinant Mouse FGF21 + Human IgG1 fusion protein (Fc Chimera)                  |
| 精製度         | > 90 % SDS-PAGE.<br>ab108580 is 0.2µm filtered                                    |
| エンドトキシン・レベル | < 0.100 Eu/µg   |
| 発現系         | HEK 293 cells   |
| タンパク質長      | Full length protein   |
| Animal free | No  |
| 由来          | Recombinant   |
| 生物種         | Mouse   |
| 予測される分子量    | 50 kDa including tags   |
| 配列の追加情報     | Mouse FGF21 (aa 1-210) is fused at the C-terminus to the Fc portion of human IgG. |
| 製品の詳細       | Recombinant Mouse FGF21 + IgG1 protein (Fc Chimera)                               |

### 特性

Our **Abpromise guarantee** covers the use of **ab108580** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

|          |  |
|----------|--|
| アプリケーション | SDS-PAGE   |
| 製品の状態    | Liquid   |
| 備考       | Working aliquots are stable for up to 3 months when stored at -20°C. |

### 前処理および保存

|            |   |
|------------|---|
| 保存方法および安定性 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.<br>Constituent: PBS |
|------------|---|

### 関連情報

|       |  |
|-------|--|
| 細胞内局在 | FGF21: Secreted. IgG1: Secreted. Cell membrane; Single-pass membrane protein |
|-------|--|

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors