

Recombinant Mouse GAL4 protein ab136705

画像数 1

製品の詳細

| | |
|-------------|--|
| 製品名 | Recombinant Mouse GAL4 protein |
| 精製度 | > 90 % SDS-PAGE. ab136705 is purified using conventional chromatography techniques. |
| 発現系 | Escherichia coli |
| アクセッション番号 | Q8K419 |
| タンパク質長 | Full length protein |
| Animal free | No |
| 由来 | Recombinant |
| 生物種 | Mouse |
| 配列 | MGSSHHHHHSSGLVPRGSHMGSMAYVPAPGYQPTYNPTLPY KRPIPGGL SVGMSVYIQMAKENMRRFHVNFVAVGQDDGADVAFHFNPRFD GWDKVVFN TMQSGQWGKEEKKKSMFPQKGGKHFELVFMVMPPEHYKVVVNGN SFYEYGHR LPVQMVTHLQVDGDLELQSFNINFLGGQPAAAPYPGAMTIPAYP AGSPGYNP PQMNTLPVMTGPPVFNPRVPYVVGALQGGLTVRRTIIKGYVL PTARNFVI NFKVGSSEGDIALHLNPRIGDSVVRNSFMNGSWGAEERKVAYN PFGPGQFF DLSIRCGMDRFKVFANGQHLFDFSHRFQAFQMVDLTLEINGDI TLSYVQI |
| 予測される分子量 | 39 kDa including tags |
| 領域 | 1 to 326 |
| タグ | His tag N-Terminus |

特性

Our **Abpromise guarantee** covers the use of **ab136705** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション SDS-PAGE

Mass Spectrometry

質量分析

MALDI-TOF

製品の状態

Liquid

前処理および保存

保存方法および安定性

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 0.02% DTT, 0.32% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

関連情報

機能

Galectin that binds lactose and a related range of sugars. May be involved in the assembly of adherens junctions.


配列類似性

Contains 2 galectin domains.

ドメイン

Contains two homologous but distinct carbohydrate-binding domains.

画像



15% SDS-PAGE analysis of 3 µg ab136705.

SDS-PAGE - Recombinant Mouse GAL4 protein
(ab136705)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors