

# Recombinant mouse Macrophage Inflammatory Protein 1 gamma ab9914

## 製品の詳細

製品名	Recombinant mouse Macrophage Inflammatory Protein 1 gamma
精製度	> 98 % SDS-PAGE. Sterile filtered Greater than 98% pure by HPLC analyses. Endotoxin level is less than 0.1 ng per g (1EU/g).
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Mouse

## 特性

Our **Abpromise guarantee** covers the use of **ab9914** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Functional Studies SDS-PAGE
製品の状態	Lyophilized
備考	The biological activity of this product is determined by its ability to chemoattract human neutrophils and human T cells using a concentration range of 1.0-10.0 ng/ml

## 前処理および保存

保存方法および安定性	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.  This product is an active protein and may elicit a biological response in vivo, handle with caution.
再構成	Please reconstitute this product using 200ul of sterile water to a concentration of 0.1mg/ml.

## 関連情報

<b>機能</b>	Monokine with inflammatory, pyrogenic and chemokinetic properties. Circulates at high concentrations in the blood of healthy animals. Binding to a high-affinity receptor activates calcium release in neutrophils. It also inhibits colony formation of bone marrow myeloid immature progenitors.
<b>組織特異性</b>	Expressed mainly in the liver, lung, and the thymus, although some expression has been detected in a wide variety of tissues except brain.
<b>配列類似性</b>	Belongs to the intercrine beta (chemokine CC) family.
<b>細胞内局在</b>	Secreted.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors