

Recombinant mouse MCP3 protein ab9912

製品の詳細

製品名	Recombinant mouse MCP3 protein
生理活性	Determined by its ability to chemoattract BaF3 cells expressing human CCR2 receptor using a concentration range of 300-500 ng/mL.
精製度	>= 98 % SDS-PAGE. = 98% by SDS-PAGE gel and HPLC analyses.s.
エンドトキシン・レベル	< 1.000 Eu/μg
発現系	Escherichia coli
アクセッション番号	<u>Q03366</u>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Mouse
配列	QPDGPNASTCCYVKKQKIPKRNLKSYRRITSSRCPWEAVIFK TKKGMEVC AEAHQKWVEEAIAYLDMKTPTPKP
予測される分子量	11 kDa
領域	24 to 97

特性

Our **Abpromise guarantee** covers the use of **ab9912** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Functional Studies
	SDS-PAGE

製品の状態	Lyophilized
-------	-------------

前処理および保存

保存方法および安定性	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. Sterile filtered through a 0.2 micron filter. Lyophilized with no additives. This product is an active protein and may elicit a biological response in vivo, handle with caution.
------------	--

再構成

For lot specific reconstitution information please contact our Scientific Support Team.

関連情報

機能	Chemotactic factor that attracts monocytes and eosinophils, but not neutrophils. Augments monocyte anti-tumor activity. Also induces the release of gelatinase B. This protein can bind heparin. Binds to CCR1, CCR2 and CCR3.
配列類似性	Belongs to the intercrine beta (chemokine CC) family.
翻訳後修飾	O-glycosylated.
細胞内局在	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors