

Recombinant Mouse PAI1 + PLAT protein ab92709

画像数 1

製品の詳細

製品名	Recombinant Mouse PAI1 + PLAT protein
精製度	> 95 % SDS-PAGE. Mouse sc-tPA is reacted with an excess of Mouse PAI1 stable mutant at physiological pH. The resulting 1:1 molar complex is purified from the remaining PAI1 by affinity chromatography with immobilized antibodies to Mouse tPA.
発現系	Insect cells
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Mouse

特性

Our **Abpromise guarantee** covers the use of **ab92709** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Liquid
備考	Solubility: > 2 mg/mL Spectrophotometric Data: $\epsilon^{0.1\%}=1.5$

前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 7.40 Constituents: 9.52% HEPES, 0.58% Sodium chloride
------------	---

関連情報


関連性	PAI1 acts as "bait" for tissue plasminogen activator (TPA), urokinase, and protein C. Its rapid interaction with TPA may function as a major control point in the regulation of fibrinolysis. TPA: This converts the abundant, but inactive, zymogen plasminogen to plasmin by hydrolyzing a single
-----	---

Arg-Val bond in plasminogen. By controlling plasmin mediated proteolysis, it plays an important role in tissue remodeling and degradation, in cell migration and many other physiopathological events.

細胞内局在

Secreted protein; extracellular space.

画像



Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors