

Recombinant mouse Prolactin/PRL protein (Active) ab54139

製品の詳細

製品名	Recombinant mouse Prolactin/PRL protein (Active)
生理活性	Determined by its ability to induce the proliferation of rat Nb2-11 cells in the concentration range of 0.1-1.0 ng/ml.
精製度	> 98 % SDS-PAGE.
発現系	Escherichia coli
アクセッション番号	P06879
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Mouse
配列	MLPICSAGDCQTSRLRFDRVVILSHYIHTLYTDMFIEFDKQ YVQDREFM VKVINDCPTSSLATPEDKEQALKVPPEVLLNLILSLVQSSSD PLFQLITG VGGIQEAPEYILSRAKEIEEQNKQLLEGVEKIISQAYPEAKG NGIYFVWS QLPSLQGVDEESKILSLRNTIRCLRRDSHKVDNFLKVLRCQI AHQNNC
予測される分子量	23 kDa
領域	30 to 226
配列の追加情報	Full length mature chain without signal peptide.

特性

Our **Abpromise guarantee** covers the use of **ab54139** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
	Functional Studies

製品の状態	Lyophilized
-------	-------------

保存方法および安定性	<p>Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.</p> <p>This product is an active protein and may elicit a biological response in vivo, handle with caution.</p>
再構成	<p>Reconstituted Prolactin is stable for at least 3 months when stored in working aliquots with a carrier protein at -20°C. Avoid repeated freeze/thaw cycles.</p>

関連情報

関連性	<p>Various hormones are secreted from the anterior pituitary during development and growth. Prolactin is a growth factor secreted by the anterior pituitary that is necessary for the proliferation and differentiation of the mammary glands.</p>
細胞内局在	<p>Secreted</p>

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors