

## Product datasheet

# Recombinant Mouse S100A4 protein ab109341

画像数 1

### 製品の詳細

製品名	Recombinant Mouse S100A4 protein
精製度	> 90 % SDS-PAGE. ab109341 is purified using conventional chromatography techniques.
発現系	Escherichia coli
アクセッション番号	<b>P07091</b>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Mouse
配列	<b>MGSSHHHHHSSGLVPRGSHMARPLEEALDIVSTFHKYS GKEGDKFKLN KTELKELLTRELPSFLGKRTDEAAFQKVMNSLDSNRDNEVDF QEYCVFLS CIAMMCNEFFEGCPDKEPRKK</b>
予測される分子量	14 kDa including tags
領域	1 to 101
タグ	His tag N-Terminus

### 特性

Our **Abpromise guarantee** covers the use of **ab109341** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE Mass Spectrometry
質量分析	MALDI-TOF
製品の状態	Liquid

### 前処理および保存

保存方法および安定性	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
------------	--

pH: 8.00

Constituents: 0.0308% DTT, 0.316% Tris HCl, 20% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

## 関連情報

---

### 組織特異性

Ubiquitously expressed.

### 配列類似性


Belongs to the S-100 family.

Contains 2 EF-hand domains.

---

## 画像

---



15% SDS-PAGE analysis of 3 µg ab109341.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors