

Recombinant Mouse SF20 protein ab54228

製品の詳細

製品名	Recombinant Mouse SF20 protein
精製度	> 95 % SDS-PAGE. Purity: > 98% by SDS-PAGE
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Mouse

特性

Our **Abpromise guarantee** covers the use of **ab54228** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Lyophilized

前処理および保存

保存方法および安定性	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
再構成	Reconstitute in water to a concentration of 0.1-1.0 mg/ml. This solution can be diluted in water or other buffer solutions.

関連情報

組織特異性	Expressed in synovial tissue. Found in synovial fluid of patients with arthropaties.
配列類似性	Belongs to the UPF0556 family.
細胞内局在	Secreted. Endoplasmic reticulum-Golgi intermediate compartment.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors