

Recombinant rabbit Urokinase protein ab92603

画像数 1

製品の詳細

製品名	Recombinant rabbit Urokinase protein	
生理活性	Activity: 200,000 IU/mg (Relative to WHO International Standard for Human High Molecular Weight Urokinase, NIBSC 87/594) Activity: 200,000 IU/mg (Relative to WHO International Standard for Human High Molecular Weight Urokinase, NIBSC 87/594)	
精製度	> 95 % SDS-PAGE.	
発現系	Insect cells	
タンパク質長	Full length protein	
Animal free	No	
由来	Recombinant	
生物種	Rabbit	
配列	MSHELHGVSDASNCGCLNGGTCVTYKYFSNIWRCNCPKKFQG EHCEIDTL KTCYHGDGHSYRGKANTDIMDRPCLAWNSANVLTKTYHAHRP DALQLGLG KHNYCRNPDHQRRPWCYVQVGLKQLIQECKVHDCSSGKKPAL PPGKLEFQ CGQKALRPRFKIIGGEFTIENQPWFAAIYRRHRGGSVTYVC GGSLISPC	
領域	21 to 433	

特性

Our Abpromise guarantee covers the use of ab92603 in the following tested applications.	
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.	
アプリケーション	Functional Studies SDS-PAGE
製品の状態	Liquid
備考	Produced using non-baculovirus insect cells.

前処理および保存

保存方法および安定性	<p>Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.</p> <p>pH: 5.00</p> <p>Constituents: 0.41% Sodium acetate, 0.0292% EDTA, 0.58% Sodium chloride</p> <p>This product is an active protein and may elicit a biological response in vivo, handle with caution.</p>
-------------------	---

関連情報

機能	Specifically cleave the zymogen plasminogen to form the active enzyme plasmin.
組織特異性	Expressed in the prostate gland and prostate cancers.
配列類似性	Belongs to the peptidase S1 family.
	Contains 1 EGF-like domain.
	Contains 1 kringle domain.
	Contains 1 peptidase S1 domain.
翻訳後修飾	Phosphorylation of Ser-158 and Ser-323 abolishes proadhesive ability but does not interfere with receptor binding.
細胞内局在	Secreted.

画像



10% SDS-PAGE showing ab92603
 Lane 1 Urokinase Non-reduced (3μg)
 Lane 2 Urokinase Reduced (3μg)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors