

Recombinant Rhesus monkey CD86 protein (His tag) ab215644

製品の詳細

製品名	Recombinant Rhesus monkey CD86 protein (His tag)
精製度	> 95 % SDS-PAGE.
エンドトキシン・レベル	< 1.000 Eu/μg
発現系	Mammalian
アクセッション番号	<u>H9ZFI8</u>
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Rhesus monkey
配列	YFNETADLPCQFANSQNRSLSELVVFQWQENLVLNEVYLGK EKFDSVHS KYMGRTSFDPESWTLRLHNLQIKDKGLYQCIHHKRPTGMIR IHQMNSEL SVLANFSQPEIVPISNITENMYINLTCSSIHGYPEPEKMSVL LRTKNSTI EYDGVMQKSQDNVTELYDVSISLSVSFPDVTSNMTIFCVLET DKTQLLSS PFSIELEDPPPPDHIPHHHHH
予測される分子量	26 kDa including tags
領域	31 to 247
タグ	His tag C-Terminus

特性

Our **Abpromise guarantee** covers the use of **ab215644** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション SDS-PAGE

製品の状態 Lyophilized

前処理および保存

保存方法および安定性 Shipped at 4°C. Store at -20°C long term. Avoid freeze / thaw cycle.

pH: 7.40
Constituent: 100% PBS

再構成 Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 µg/ml. Dissolve the lyophilized protein in ddH₂O.

関連情報

機能	Receptor involved in the costimulatory signal essential for T-lymphocyte proliferation and interleukin-2 production, by binding CD28 or CTLA-4. May play a critical role in the early events of T-cell activation and costimulation of naive T-cells, such as deciding between immunity and anergy that is made by T-cells within 24 hours after activation. Isoform 2 interferes with the formation of CD86 clusters, and thus acts as a negative regulator of T-cell activation.
組織特異性	Expressed by activated B-lymphocytes and monocytes.
配列類似性	Contains 1 Ig-like C2-type (immunoglobulin-like) domain. Contains 1 Ig-like V-type (immunoglobulin-like) domain.
翻訳後修飾	Polyubiquitinated; which is promoted by MARCH8 and results in endocytosis and lysosomal degradation.
細胞内局在	Membrane.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors