

# Recombinant Toxoplasma gondii GRA 1 protein ab43137

## 製品の詳細

製品名	Recombinant Toxoplasma gondii GRA 1 protein
生理活性	Immunoreactive with sera of Toxoplasma gondii-infected individuals.
精製度	> 95 % SDS-PAGE. Sephacrose-Derived Purification determined by 10% PAGE (coomassie staining).
発現系	Escherichia coli
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
タグ	GST tag N-Terminus
配列の追加情報	The protein contains the P24 (GRA1) immunodominant region.

## 特性

Our **Abpromise guarantee** covers the use of **ab43137** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	ELISA Western blot
製品の状態	Liquid
備考	E.coli derived recombinant protein that contain p24 (GRA1) region, uniprot P13403.

## 前処理および保存

保存方法および安定性	Shipped at 4°C. Store at -20°C. Constituents: 9% Urea, 0.79% Tris HCl, 50% Glycerol
------------	--

## 関連情報

関連性	GRA 1 is a major antigen secreted by Toxoplasma gondii.
細胞内局在	Secreted; Located in dense granules of tachyzoites.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors