

# Recombinant VZV ORF26 protein ab48958

### 製品の詳細

---

製品名	Recombinant VZV ORF26 protein
精製度	> 90 % SDS-PAGE.
発現系	Escherichia coli
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant

### 特性

---

Our **Abpromise guarantee** covers the use of **ab48958** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	ELISA
	Western blot

製品の状態	Liquid
-------	--------

### 前処理および保存

---

保存方法および安定性	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. pH: 7.2 Preservative: 0.01% Sodium azide Constituents: 9% Urea, 0.395% Tris HCl, 0.0292% EDTA, 50% Glycerol (glycerin, glycerine)
------------	--

### 関連情報

---

関連性	Varicella-zoster virus (VZV or HHV3) is a member of the genus Varicellovirus in the Alphaherpesvirinae subfamily of the Herpesviridae. It is the causative agent of chicken pox (varicella) in children, after which it establishes latency in the sensory ganglia with the potential to reactivate at a later time to cause shingles (zoster). This is an extremely stable virus. The genome is comprised of ~125 kb of linear double-stranded DNA containing approximately 71 open reading frames (ORFs). The viral structure is similar to that of other alphaherpesviruses, consisting of two unique regions, unique long and unique short, each flanked by inverted repeats; short repeats termed terminal repeat long and internal repeat long border the unique long region,
-----	---

while larger repeats termed terminal repeat short (TRS) and internal repeat short (IRS) border the unique short region. Varicella-zoster virus (VZV) interacts with cell surface heparan sulfate proteoglycans during virus attachment. VZV-gE is a glycoprotein that plays an active or supportive role in VZV cell membrane fusion. VZV-gE was found to enhance the fusogenic potential of VZV gB.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors