

Sample Diluent NS ab193972

SimpleStep ELISA[®]

1 References

製品の概要

製品名	Sample Diluent NS
製品の概要	<p>Abcam's Sample Diluent NS is a proprietary reagent diluting samples to be used in immunoassays.</p> <p>This diluent is specially formulated to quantitate your proteins of interest in complex sample types, such as Human serum or Human plasma.</p>
特記事項	<p><u>Instructions for use</u></p> <p>Sample Diluent NS comes as a ready to use solution.</p> <p>Samples should be diluted to within the desired working range of your assay in Sample Diluent NS.</p>

製品の特性

保存方法 Store at +4°C. Please refer to protocols.

内容	12 ml	50 ml	100 ml
Sample Diluent NS (ab193972)	1 x 12ml	1 x 50ml	1 x 100ml

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors