

Chicken Anti-Goat IgG H&L ab86245

[2 References](#) [画像数 1](#)

製品の概要

製品名	Chicken Anti-Goat IgG H&L
由来種	Chicken
ターゲット生物種	Goat
特異性	The antibody reacts with both heavy and light chains of goat IgG, and may have some cross-reactivity with IgG of other species.
アプリケーション	適用あり: WB, ELISA
免疫原	Purified Goat IgG.
ポジティブ・コントロール	Goat plasma.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.02% Sodium azide Constituent: PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgY

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab86245の使用に適用されず**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		
ELISA		

追加情報

ELISA: 1/2000 - 1/10000.

WB: 1/2000 - 1/10000 for color substrates.

WB: 1/10000 - 1/40000 for chemiluminescence.

Detects a band of approximately 30 and 60 kDa.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

画像



Chicken Anti-Goat IgG H&L (ab86245) at 1/4000 dilution + Goat plasma at 5 µg

Observed band size : 30 + 60 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors