

Chicken Anti-Rat IgG H&L (HRP) ab6836

6 References

製品の概要

| | |
|----------|---|
| 製品名 | Chicken Anti-Rat IgG H&L (HRP) |
| 由来種 | Chicken |
| ターゲット生物種 | Rat |
| アプリケーション | 適用あり: WB, ELISA, Immunomicroscopy, Dot blot, ICC/IF, IHC-P, IHC-Fr |
| 免疫原 | Rat IgG whole molecule |
| 標識 | HRP |

製品の特性

| | |
|----------|---|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Store at +4°C. |
| バッファー | Preservative: 0.01% Gentamicin sulphate Constituents: 0.878% Sodium chloride, 1% BSA, 0.424% Potassium phosphate |
| 精製度 | Affinity purified |
| 特記事項(精製) | Horseradish Peroxidase (HRP). |
| 特記事項(標識) | This product was prepared from monospecific antiserum by immunoaffinity chromatography using Rat IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities |
| ポリ/モノ | ポリクローナル |

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab6836の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|-------------------------------------|
| WB | | Use at an assay dependent dilution. |
| ELISA | | Use at an assay dependent dilution. |

| アプリケーション | Abreviews | 特記事項 |
|------------------|-----------|--|
| Immunomicroscopy | | Use at an assay dependent concentration. |
| Dot blot | | Use at an assay dependent concentration. |
| ICC/IF | | 1/500 - 1/2500. |
| IHC-P | | Use at an assay dependent concentration. |
| IHC-Fr | | Use at an assay dependent concentration. |

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors