

# Donkey Anti-Goat IgG H&L (Alkaline Phosphatase) preadsorbed ab7126

### 製品の概要

|          |  |
|----------|--|
| 製品名      | Donkey Anti-Goat IgG H&L (Alkaline Phosphatase) preadsorbed                          |
| 由来種      | Donkey   |
| ターゲット生物種 | Goat   |
| アプリケーション | <b>適用あり:</b> ICC/IF, Dot blot, ELISA, IHC-P, IHC-Fr, Immunomicroscopy, WB            |
| 吸着処理血清   | Chicken, Guinea pig, Hamster, Horse, Mouse, Rabbit, Rat <a href="#">more details</a> |
| 免疫原      | Goat IgG whole molecule  |
| 標識       | Alkaline Phosphatase   |

### 製品の特性

|           |  |
|-----------|--|
| 製品の状態     | Liquid   |
| 保存方法      | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.                                 |
| バッファー     | Preservative: 0.01% Sodium azide<br>Constituents: 0.88% Sodium chloride, 1% BSA, 50% Glycerol, 0.0014% Zinc chloride, 0.0095% Magnesium chloride, 0.79% Tris HCl |
| 精製度       | Affinity purified  |
| 特記事項 (精製) | This product was prepared from monospecific antiserum by immunoaffinity chromatography using Goat IgG coupled to agarose beads.                                  |
| 特記事項 (標識) | Alkaline Phosphatase (Calf Intestine) (Molecular Weight 140,000 daltons)   |
| ポリ/モノ     | ポリクローナル  |

### アプリケーション

**The Abpromise guarantee** Abpromise保証は、次のテスト済みアプリケーションにおけるab7126の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

| アプリケーション         | Abreviews | 特記事項                                     |
|------------------|-----------|--|
| ICC/IF           |           | 1/200 - 1/1000.                          |
| Dot blot         |           | Use at an assay dependent concentration. |
| ELISA            |           | 1/3000 - 1/12000.                        |
| IHC-P            |           | Use at an assay dependent concentration. |
| IHC-Fr           |           | Use at an assay dependent concentration. |
| Immunomicroscopy |           | Use at an assay dependent concentration. |
| WB               |           | 1/500 - 1/2000.                          |

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors