

Donkey Anti-Goat IgG H&L (DyLight® 550) preadsorbed ab96936

3 References [画像数 1](#)

製品の概要

製品名	Donkey Anti-Goat IgG H&L (DyLight® 550) preadsorbed
由来種	Donkey
ターゲット生物種	Goat
特異性	By immunoelectrophoresis and ELISA this antibody reacts specifically with Goat IgG and with light chains common to other Goat immunoglobulins. No antibody was detected against non immunoglobulin serum proteins. Reduced cross-reactivity to chicken, horse, human, mouse, pig, rabbit and rat was detected.
アプリケーション	適用あり: IHC-P, ICC/IF, Flow Cyt
吸着処理血清	Chicken, Horse, Human, Mouse, Pig, Rabbit, Rat more details
免疫原	Other Immunogen Type corresponding to Goat IgG.
標識	DyLight® 550. Ex: 562nm, Em: 576nm

製品の特性


製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: 0.2% BSA, PBS
精製度	Immunogen affinity purified
特記事項(精製)	Antiserum was cross adsorbed using chicken, human, mouse, pig, rabbit and rat immunosorbents to remove cross reactive antibodies. This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to DyLight® 550.
ポリモノ	ポリクローナル
アイソタイプ	IgG
特記事項	DyLight® 550 replaces DyLight® 549.

アプリケーション

The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab96936の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P		1/50 - 1/500. DyLight® 550 replaces DyLight® 549. The DyLight® 550 excitation is 562nm and the emission is 576nm. DyLight® 550 provides the same orange-to-red fluorescence and photostability as the DyLight® 549 dye. Please refer to the vial in order to check whether the product you received is DyLight® 549 or DyLight® 550. This information however will have no impact on your
ICC/IF		1/50 - 1/500.
Flow Cyt		1/50 - 1/200.

画像



Emission spectra of DyLight® fluorochromes available in our catalog.

Line colors represent the approximate visible colors of the wavelength maxima.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors