

Donkey Anti-Mouse IgG H&L (DyLight® 488) ab96875

18 References [画像数 2](#)

製品の概要

製品名	Donkey Anti-Mouse IgG H&L (DyLight® 488)
由来種	Donkey
ターゲット生物種	Mouse
特異性	By immunoelectrophoresis and ELISA this antibody reacts specifically with Mouse IgG and with light chains common to other Mouse immunoglobulins. No antibody was detected against non immunoglobulin serum proteins.
アプリケーション	適用あり: ICC/IF, IHC-P, Flow Cyt
免疫原	Other Immunogen Type corresponding to Mouse IgG.
標識	DyLight® 488. Ex: 493nm, Em: 518nm

製品の特性


製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 6.8 Preservative: 0.09% Sodium azide Constituents: 0.2% BSA, PBS
精製度	Immunogen affinity purified
特記事項 (精製)	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to DyLight® 488.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab96875の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC/IF		1/50 - 1/500.
IHC-P		1/50 - 1/500.
Flow Cyt		1/50 - 1/200.

画像




Immunocytochemistry/ Immunofluorescence -
Donkey Anti-Mouse IgG H&L (DyLight® 488)
(ab96875)

This image is courtesy of an Abreview submitted by Charles Pallangyo.

ab7817 staining alpha smooth muscle actin in mouse primary colon myofibroblasts by ICC/IF

(Immunocytochemistry/immunofluorescence). Cells were fixed with acetone and blocked with 5% BSA for 30 minutes at 25°C. Samples were incubated with primary antibody (1/100 in PBS + 5% BSA) for 2 hour at 25°C. A DyLight® 488-conjugated donkey anti-mouse IgG H&L (ab96875) (1/1000) was used as the secondary antibody. Vimentin appears in red.



Donkey Anti-Mouse IgG H&L (DyLight® 488)
(ab96875)

Emission spectra of DyLight® fluorochromes available in our catalog.

Line colors represent the approximate visible colors of the wavelength maxima.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors